



Controller 2.0

Smart
Reliable
Accurate

Our next-generation controller delivers class-leading ease of use, giving your patients greater freedom with fewer battery changes and associated alarms. In addition, it provides easily accessible and actionable data to help you streamline patient management.

Controller 2.0 is the latest advancement demonstrating HeartWare's commitment to continuous improvement.

HeartWare®

What's new in Controller 2.0?

Smarter use of battery power

Our new battery drain scheme allows patients to go up to 14 hours before performing a battery exchange, while still providing generous time to exchange power sources.

Easier patient management

Now you can easily view battery cycle counts, waveform amplitude, and pump settings directly from Controller 2.0—without the need for a monitor.

Upgraded connections

Controller 2.0's new metal connectors are designed to increase wear resistance. They also now provide a natural, 12 o'clock orientation to give you and your patients more reliable connections.

Ability to begin recording data...from the beginning

The system now begins recording logfile data less than a minute after powering on to support your team's investigational work.

Updated controller, familiar look

Over 50% of Controller 2.0 components have been replaced and over 30 electrical, mechanical, and software reliability-focused updates were made. The combination of new components and dozens of software changes are designed to make the HeartWare HVAD® System as robust as ever.

Compatible with older electrical systems and generators

Controller 2.0 no longer requires a 3-pronged, grounded outlet for the AC adapter, allowing it to be compatible with older electrical systems and generators.

14 Controller languages to meet your patients' needs

You can now customize Controller 2.0 to view data in 14 different languages, allowing you to provide easy-to-understand guidance to even more patients.

Available languages: English, Danish, Dutch, Finnish, French, German, Italian, Japanese*, Norwegian, Polish, Portuguese*, Spanish, Swedish and Turkish*

*New to International

HeartWare

HeartWare, Inc.
14400 NW 60th Avenue
Miami Lakes, FL 33014 USA
www.heartware.com

HEARTWARE, HVAD, and the HEARTWARE logo are trademarks of HeartWare, Inc.

© 2016 HeartWare, Inc. INTL1080 Rev01 08/16

CE